



Automatic Approach Seeking Optimal Frequency Modulation Parameters In Chirp Inversion and Chirp Reversal Ultrasound Contrast Imaging

Amira Zaylaa, Sébastien Ménigot, Jamal Charara, Jean-Marc Girault

► To cite this version:

Amira Zaylaa, Sébastien Ménigot, Jamal Charara, Jean-Marc Girault. Automatic Approach Seeking Optimal Frequency Modulation Parameters In Chirp Inversion and Chirp Reversal Ultrasound Contrast Imaging. The Second Science Meeting. Doctoral School of Sciences and Technology, Lebanese University, (2012)., Jun 2012, Hadath, Lebanon. inserm-01182641

HAL Id: inserm-01182641

<https://www.hal.inserm.fr/inserm-01182641>

Submitted on 1 Aug 2015

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Automatic Approach Seeking Optimal Frequency Modulation Parameters In Chirp Inversion and Chirp Reversal Ultrasound Contrast Imaging

A. Zaylaa^{a,b,c}, S. Ménigot^{a,b}, J. Charara^d and J.-M. Girault^{a,b}

a. Imaging and Brain Laboratory, Université François Rabelais de Tours, B.P. 22, France. amira.zaylaa@etu.univ-tours.fr

b. INSERM U930, Tours, France.

c. P.R.A.S.E., Doctoral School of Sciences and Technology, Lebanese University, Hadath, Lebanon.

d. Departement of Physics and Electronics, Lebanese University, Hadath, Lebanon.

I – Introduction


State-of-the art applications involving frequency coded excitations known as chirp excitations, are cutting edge in contrast ultrasound imaging. A question may arise: How does the manufacturer select the frequency modulation parameters of chirp excitations in clinical practice? The present study aims to select, automatically, the optimal frequency modulation parameters of chirp excitations by maximizing the backscattered power using two post-processing techniques: chirp inversion and chirp reversal imaging. Why the power? As the contrast-to-tissue ratio \nearrow the power backscattered by microbubbles \nearrow

$$CTR = \frac{P_{microbubbles}}{P_{tissues}}$$

II – Materials and Methods

Simulated Microbubbles with Matlab

- Air-filled of diameter 5 μm
- Immersed in blood mimicking fluid
- Motion described by Modified Plesset-Rayleigh Differential Equation




Emission/Reception

- Transducer: Bandwidth 63%


Chirp Excitations

- Frequency modulated \rightarrow Linear frequency modulation (FM) law
- Initial driving pressure \rightarrow 100 kPa \rightarrow NonLinear Response
- Envelope \rightarrow Gaussian


Transmitted Chirp
($\beta=0.03\text{MHz/s}$, $f_0=1.5\text{ MHz}$)




Chirp Spectrogram




Block Diagram of Chirp Reversal and/or Chirp Inversion post-processing techniques.



Chirp Reversal



Chirp Inversion



Automatic Optimization Approach

- Cost Function J \rightarrow Resultant Backscattered Power


$$w_{k+1} = w_k + \mu_k \cdot \nabla J(w_k)$$

- Automatic Optimization: Maximum Number of iterations k=50

$$\begin{pmatrix} f_{k+1} \\ \beta_{k+1} \end{pmatrix} = \begin{pmatrix} f_k \\ \beta_k \end{pmatrix} + \begin{pmatrix} \mu_{11} & \mu_{12} \\ \mu_{21} & \mu_{22} \end{pmatrix} \times \begin{pmatrix} \nabla f^P \\ \nabla \beta^P \end{pmatrix} \Rightarrow \text{Optimum}(f_o, \beta)$$

III – Simulation Results

- Automatic versus Standard chirp reversal (currently available in literature) Gain \rightarrow 5.4 dB
- Automatic versus Standard chirp inversion Gain \rightarrow 2.6 dB
- Resultant backscattered power attained after 15 iterations
- Central Transmitted frequency \nearrow 1.1 MHz
- Frequency modulation index \nearrow 0.23 MHz/s



Processing Technique	Chirp Reversal		Chirp Inversion	
	Optimization	Stand.	Auto.	Stand.
Modulation Law	Linear FM	Linear FM	Linear FM	Linear FM
Gain (dB)	-24.3	-41.7	-24.5	-31.5
f_0 (MHz)	2.2	3	2.2	2.75
Beta (MHz/s)	0.23	0.2	0.23	0.1

IV – Conclusion

- Automatic Optimization by chirp reversal/ inversion converged to similar optimized parameters
- Power gain by Chirp Reversal was advantageous over chirp inversion
- Power Gain attained by automatic chirp reversal was significant compared to standard technique
- The proposed Automatic approach may improve CTR of Ultrasound Images so we tend to validate it by experiment